
OSOM[®] Strep A Test 141E

NO

OSOM® Strep A Test

Katalognummer 141E eller 141E-20

CLIA-kompleksitet: godkjent

TILTENKT BRUK

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No. 141E - OSOM Strep A Test (50 Tests)

No. 141E-20 - OSOM Strep A Test (20 Tests)

1. Youmans, G.P., Paterson, P.Y. and Sommers, H.M. Upper Respiratory Tract Infection: General Considerations, in The Biologic and Clinical Basis of Infectious Diseases, W. B. Saunders Co., Philadelphia, 177-183, 1980.
2. Faklam, R.R. and Washington, J.A., Streptococcus and Related Catalase-Negative Gram-Positive Cocci, in Manual of Clinical Microbiology, 5th Edition, Balows, A., Hausler, W.J., Herrmann, K.L., Isenberg, H.D., and Shadomy, H.J., Eds., Am. Society of Microbiology, Washington, D.C., 238-257, 1991.
3. CDC, Biosafety in Microbiological and Biomedical Laboratories, 2nd Ed., HHS Publication No. 8808395, 4-6, 1988.
4. American Academy of Pediatrics. Summaries of Infectious Disease. In: Pickering LK, ed. 2000 Red Book: Report of the Committee on Infectious Diseases. 25th ed. Elk Grove Village IL: American Academy of Pediatrics; 2000 528
5. Lauer, B.A., Reller, L.B., and Mirrett, S., Effect of Atmosphere and Duration of Incubation on Primary Isolation of Group A Streptococci from Throat Cultures, J. Clin. Microb., 17: 338-340, 1983.
6. Wannamaker, L.W., Differences Between Streptococcal Infections of the Throat and of the Skin, N. Eng. J. Med., 282: 23-31, 78-85, 1970

DEFINITIONS OF SYMBOLS / SYMBOLFORKLARING / DÉFINITIONS DES SYMBOLES / DEFINIZIONE DEI SIMBOLI / DEFINITIONEN DER SYMBOLE / SYMBOLER / DEFINICIONES DE SÍMBOLOS / SYMBOLDEFINITIONER / DEFINIÇÕES DOS SÍMBOLOS / OBJAŚNIENIE SYMBOLI

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
**In Vitro Diagnostic Medical Device**

Medicinsk udstyr til in vitro-diagnostik
Dispositif médical de diagnostic in vitro
Dispositivo medico-diagnostico in vitro
In-vitro-Diagnostikum
Til in vitro diagnostisk bruk
Producto sanitario para diagnóstico in vitro
Medicintekniska produkter för in vitro diagnostik
Dispositivo médico para diagnóstico in vitro
Wyrób do diagnostyki In Vitro

**Toxic**

Giftig
Toxique
Tossico
Giftig
Giftig
Tóxico
Giftig
Tóxico
Toksyczne

**Authorised representative in the EC**

Repræsentant i det Europæiske Fællesskab
Mandataire dans la Communauté européenne
Mandatario nella Comunità Europea
Bevollmächtigter in der Europäischen Gemeinschaft
Autorisert representant
Representante autorizado en la Comunidad Europea
Auktoriserad representant i Europeiska Gemenskapen
Representante autorizado na comunidade Europeia
Autoryzowany przedstawiciel w Unii Europejskiej

**Contains sufficient for <n> tests**

Indeholder tilstrækkeligt til "n" test
Contenu suffisant pour "n" tests
Contenuto sufficiente per "n" saggi
Inhalt ausreichend für <n> Prüfungen
Antall tester
Contenido suficiente para <n> ensayos
Räcker till "n" antal tester
Conteúdo suficiente para <n> ensaios
Zawiera ilość wystarczającą do wykonania <n> badań

OSOM and the QC Inside logo

are registered U.S. trademarks of Sekisui Diagnostics, LLC.

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Licensed under U.S. Patent Nos. 5,714,389; 5,989,921; 6,485,982 and 6,979,576 and related non-U.S. patents and patent applications.

Autoriseret under U.S.A. Patentnumre 5,714,389, 5,989,921, 6,485,982 og 6,979,576 samt relaterede ikke-U.S.A. patenter og patentansøgninger.

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Sekisui Diagnostics, LLC
6659 Top Gun Street
San Diego, CA 92121
USA

Sekisui Diagnostics (UK) Limited
50 Gibson Drive
Kings Hill, West Malling
Kent, ME19 4AF, UK

SEKISUI
DIAGNOSTICS

Framingham, MA 01701
Tel: 800-332-1042
www.sekisuidiagnostics.com

© 2011 Sekisui Diagnostics, LLC – All rights reserved.

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sind eingetragene US-Marken von Sekisui Diagnostics, LLC
er registrerte varemerker i USA for Sekisui Diagnostics, LLC
son marcas registradas de EE.UU. Sekisui Diagnostics, LLC
är registrerade varumärken i USA av Sekisui Diagnostics, LLC
são registrados as marcas de comércio de Estados Unidos do Sekisui Diagnostics, LLC
są zarejestrowanymi znakami towarowymi fi rmy Sekisui Diagnostics, LLC